

Well #	Salt	Well #	Buffer ◇
1. (A1)	1.8 M Sodium acetate trihydrate pH 7.0	1. (A1)	0.1 M BIS-TRIS propane pH 7.0
2. (A2)	2.8 M Sodium acetate trihydrate pH 7.0	2. (A2)	0.1 M BIS-TRIS propane pH 7.0
3. (A3)	1.5 M Ammonium chloride	3. (A3)	0.1 M Sodium acetate trihydrate pH 4.6
4. (A4)	1.5 M Ammonium chloride	4. (A4)	0.1 M BIS-TRIS propane pH 7.0
5. (A5)	1.5 M Ammonium chloride	5. (A5)	0.1 M Tris pH 8.5
6. (A6)	3.5 M Ammonium chloride	6. (A6)	0.1 M Sodium acetate trihydrate pH 4.6
7. (A7)	3.5 M Ammonium chloride	7. (A7)	0.1 M BIS-TRIS propane pH 7.0
8. (A8)	3.5 M Ammonium chloride	8. (A8)	0.1 M Tris pH 8.5
9. (A9)	2.2 M Sodium chloride	9. (A9)	0.1 M Sodium acetate trihydrate pH 4.6
10. (A10)	2.2 M Sodium chloride	10. (A10)	0.1 M BIS-TRIS propane pH 7.0
11. (A11)	2.2 M Sodium chloride	11. (A11)	0.1 M Tris pH 8.5
12. (A12)	3.2 M Sodium chloride	12. (A12)	0.1 M Sodium acetate trihydrate pH 4.6
13. (B1)	3.2 M Sodium chloride	13. (B1)	0.1 M BIS-TRIS propane pH 7.0
14. (B2)	3.2 M Sodium chloride	14. (B2)	0.1 M Tris pH 8.5
15. (B3)	1.0 M Ammonium citrate dibasic	15. (B3)	0.1 M Sodium acetate trihydrate pH 4.6
16. (B4)	1.8 M Ammonium citrate dibasic	16. (B4)	0.1 M Sodium acetate trihydrate pH 4.6
17. (B5)	1.0 M Ammonium citrate tribasic pH 7.0	17. (B5)	0.1 M BIS-TRIS propane pH 7.0
18. (B6)	2.0 M Ammonium citrate tribasic pH 7.0	18. (B6)	0.1 M BIS-TRIS propane pH 7.0
19. (B7)	0.7 M Sodium citrate tribasic dihydrate	19. (B7)	0.1 M BIS-TRIS propane pH 7.0
20. (B8)	0.7 M Sodium citrate tribasic dihydrate	20. (B8)	0.1 M Tris pH 8.5
21. (B9)	1.2 M Sodium citrate tribasic dihydrate	21. (B9)	0.1 M BIS-TRIS propane pH 7.0
22. (B10)	1.2 M Sodium citrate tribasic dihydrate	22. (B10)	0.1 M Tris pH 8.5
23. (B11)	0.4 M Magnesium formate dihydrate	23. (B11)	0.1 M Sodium acetate trihydrate pH 4.6
24. (B12)	0.4 M Magnesium formate dihydrate	24. (B12)	0.1 M BIS-TRIS propane pH 7.0
25. (C1)	0.4 M Magnesium formate dihydrate	25. (C1)	0.1 M Tris pH 8.5
26. (C2)	0.7 M Magnesium formate dihydrate	26. (C2)	0.1 M BIS-TRIS propane pH 7.0
27. (C3)	2.0 M Sodium formate	27. (C3)	0.1 M Sodium acetate trihydrate pH 4.6
28. (C4)	2.0 M Sodium formate	28. (C4)	0.1 M BIS-TRIS propane pH 7.0
29. (C5)	2.0 M Sodium formate	29. (C5)	0.1 M Tris pH 8.5
30. (C6)	3.5 M Sodium formate	30. (C6)	0.1 M Sodium acetate trihydrate pH 4.6
31. (C7)	3.5 M Sodium formate	31. (C7)	0.1 M BIS-TRIS propane pH 7.0
32. (C8)	3.5 M Sodium formate	32. (C8)	0.1 M Tris pH 8.5
33. (C9)	1.2 M DL-Malic acid pH 7.0	33. (C9)	0.1 M BIS-TRIS propane pH 7.0
34. (C10)	2.2 M DL-Malic acid pH 7.0	34. (C10)	0.1 M BIS-TRIS propane pH 7.0
35. (C11)	1.4 M Sodium malonate pH 7.0	35. (C11)	0.1 M BIS-TRIS propane pH 7.0
36. (C12)	2.4 M Sodium malonate pH 7.0	36. (C12)	0.1 M BIS-TRIS propane pH 7.0
37. (D1)	2.5 M Ammonium nitrate	37. (D1)	0.1 M Sodium acetate trihydrate pH 4.6
38. (D2)	2.5 M Ammonium nitrate	38. (D2)	0.1 M BIS-TRIS propane pH 7.0
39. (D3)	2.5 M Ammonium nitrate	39. (D3)	0.1 M Tris pH 8.5
40. (D4)	6.0 M Ammonium nitrate	40. (D4)	0.1 M Sodium acetate trihydrate pH 4.6
41. (D5)	6.0 M Ammonium nitrate	41. (D5)	0.1 M BIS-TRIS propane pH 7.0
42. (D6)	6.0 M Ammonium nitrate	42. (D6)	0.1 M Tris pH 8.5
43. (D7)	1.5 M Sodium nitrate	43. (D7)	0.1 M Sodium acetate trihydrate pH 4.6
44. (D8)	1.5 M Sodium nitrate	44. (D8)	0.1 M BIS-TRIS propane pH 7.0
45. (D9)	1.5 M Sodium nitrate	45. (D9)	0.1 M Tris pH 8.5
46. (D10)	4.0 M Sodium nitrate	46. (D10)	0.1 M Sodium acetate trihydrate pH 4.6
47. (D11)	4.0 M Sodium nitrate	47. (D11)	0.1 M BIS-TRIS propane pH 7.0
48. (D12)	4.0 M Sodium nitrate	48. (D12)	0.1 M Tris pH 8.5

◇ Buffer pH is that of a 1.0 M stock prior to dilution with other reagent components: pH with HCl or NaOH.

*SaltRx HT (Deep Well Block) contains ninety-six unique reagents beginning at position A1. To determine the formulation of each reagent, simply read across the page.*

Well #	Salt	Well #	Buffer ◇
49.(E1)	1.0 M Ammonium phosphate monobasic	49.(E1)	0.1 M Sodium acetate trihydrate pH 4.6
50.(E2)	1.8 M Ammonium phosphate monobasic	50.(E2)	0.1 M Sodium acetate trihydrate pH 4.6
51.(E3)	1.5 M Ammonium phosphate dibasic	51.(E3)	0.1 M Tris pH 8.5
52.(E4)	2.4 M Ammonium phosphate dibasic	52.(E4)	0.1 M Tris pH 8.5
53.(E5)	1.0 M Sodium phosphate monobasic monohydrate, Potassium phosphate dibasic / pH 5.0	53.(E5)	None
54.(E6)	1.0 M Sodium phosphate monobasic monohydrate, Potassium phosphate dibasic / pH 6.9	54.(E6)	None
55.(E7)	1.0 M Sodium phosphate monobasic monohydrate, Potassium phosphate dibasic / pH 8.2	55.(E7)	None
56.(E8)	1.8 M Sodium phosphate monobasic monohydrate, Potassium phosphate dibasic / pH 5.0	56.(E8)	None
57.(E9)	1.8 M Sodium phosphate monobasic monohydrate, Potassium phosphate dibasic / pH 6.9	57.(E9)	None
58.(E10)	1.8 M Sodium phosphate monobasic monohydrate, Potassium phosphate dibasic / pH 8.2	58.(E10)	None
59.(E11)	0.5 M Succinic acid pH 7.0	59.(E11)	0.1 M BIS-TRIS propane pH 7.0
60.(E12)	1.0 M Succinic acid pH 7.0	60.(E12)	0.1 M BIS-TRIS propane pH 7.0
61.(F1)	1.5 M Ammonium sulfate	61.(F1)	0.1 M Sodium acetate trihydrate pH 4.6
62.(F2)	1.5 M Ammonium sulfate	62.(F2)	0.1 M BIS-TRIS propane pH 7.0
63.(F3)	1.5 M Ammonium sulfate	63.(F3)	0.1 M Tris pH 8.5
64.(F4)	2.5 M Ammonium sulfate	64.(F4)	0.1 M Sodium acetate trihydrate pH 4.6
65.(F5)	2.5 M Ammonium sulfate	65.(F5)	0.1 M BIS-TRIS propane pH 7.0
66.(F6)	2.5 M Ammonium sulfate	66.(F6)	0.1 M Tris pH 8.5
67.(F7)	0.8 M Lithium sulfate monohydrate	67.(F7)	0.1 M Sodium acetate trihydrate pH 4.6
68.(F8)	0.8 M Lithium sulfate monohydrate	68.(F8)	0.1 M BIS-TRIS propane pH 7.0
69.(F9)	0.8 M Lithium sulfate monohydrate	69.(F9)	0.1 M Tris pH 8.5
70.(F10)	1.5 M Lithium sulfate monohydrate	70.(F10)	0.1 M Sodium acetate trihydrate pH 4.6
71.(F11)	1.5 M Lithium sulfate monohydrate	71.(F11)	0.1 M BIS-TRIS propane pH 7.0
72.(F12)	1.5 M Lithium sulfate monohydrate	72.(F12)	0.1 M Tris pH 8.5
73.(G1)	1.0 M Magnesium sulfate hydrate	73.(G1)	0.1 M Sodium acetate trihydrate pH 4.6
74.(G2)	1.0 M Magnesium sulfate hydrate	74.(G2)	0.1 M BIS-TRIS propane pH 7.0
75.(G3)	1.0 M Magnesium sulfate hydrate	75.(G3)	0.1 M Tris pH 8.5
76.(G4)	1.8 M Magnesium sulfate hydrate	76.(G4)	0.1 M Sodium acetate trihydrate pH 4.6
77.(G5)	1.8 M Magnesium sulfate hydrate	77.(G5)	0.1 M BIS-TRIS propane pH 7.0
78.(G6)	1.8 M Magnesium sulfate hydrate	78.(G6)	0.1 M Tris pH 8.5
79.(G7)	0.7 M Ammonium tartrate dibasic	79.(G7)	0.1 M Sodium acetate trihydrate pH 4.6
80.(G8)	0.7 M Ammonium tartrate dibasic	80.(G8)	0.1 M BIS-TRIS propane pH 7.0
81.(G9)	0.7 M Ammonium tartrate dibasic	81.(G9)	0.1 M Tris pH 8.5
82.(G10)	1.0 M Ammonium tartrate dibasic	82.(G10)	0.1 M Sodium acetate trihydrate pH 4.6
83.(G11)	1.3 M Ammonium tartrate dibasic	83.(G11)	0.1 M BIS-TRIS propane pH 7.0
84.(G12)	1.4 M Ammonium tartrate dibasic	84.(G12)	0.1 M Tris pH 8.5
85.(H1)	0.6 M Potassium sodium tartrate tetrahydrate	85.(H1)	0.1 M BIS-TRIS propane pH 7.0
86.(H2)	1.2 M Potassium sodium tartrate tetrahydrate	86.(H2)	0.1 M BIS-TRIS propane pH 7.0
87.(H3)	0.6 M Potassium sodium tartrate tetrahydrate	87.(H3)	0.1 M Tris pH 8.5
88.(H4)	1.2 M Potassium sodium tartrate tetrahydrate	88.(H4)	0.1 M Tris pH 8.5
89.(H5)	0.5 M Potassium thiocyanate	89.(H5)	0.1 M Sodium acetate trihydrate pH 4.6
90.(H6)	0.5 M Potassium thiocyanate	90.(H6)	0.1 M BIS-TRIS propane pH 7.0
91.(H7)	0.5 M Potassium thiocyanate	91.(H7)	0.1 M Tris pH 8.5
92.(H8)	4.0 M Ammonium acetate	92.(H8)	0.1 M Sodium acetate trihydrate pH 4.6
93.(H9)	4.0 M Ammonium acetate	93.(H9)	0.1 M BIS-TRIS propane pH 7.0
94.(H10)	4.0 M Ammonium acetate	94.(H10)	0.1 M Tris pH 8.5
95.(H11)	35% v/v Tacsimate pH 7.0	95.(H11)	0.1 M BIS-TRIS propane pH 7.0
96.(H12)	60% v/v Tacsimate pH 7.0	96.(H12)	0.1 M BIS-TRIS propane pH 7.0

◇ Buffer pH is that of a 1.0 M stock prior to dilution  
with other reagent components: pH with HCl or NaOH.

*SaltRx HT (Deep Well Block) contains ninety-six unique reagents beginning at position A1.  
To determine the formulation of each reagent, simply read across the page.*

34 Journey

Aliso Viejo, CA 92656-3317 U.S.A.

Tel: (949) 425-1321 • Fax: (949) 425-1611

E-mail: tech@hrmail.com

Website: www.hamptonresearch.com

**HAMPTON**  
RESEARCH

*Solutions for Crystal Growth*

© 1991-2011 Hampton Research Corp. all rights reserved  
Printed in the United States of America. This guide or parts thereof may not be reproduced in any form without the written permission of the publishers.